

Product Datasheet

GRM6 peptide (orb217370)

Description	The peptide is used to block anti-mGLUR6 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GRM6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214014.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35349 , Q5NCH9 , O15303
Entrez	108072 , 2916 , 24419
Expiration Date	12 months from date of receipt.