

Product Datasheet

GRM7 peptide (orb217371)

Description	The peptide is used to block anti-mGLUR7 antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	GRM7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214015.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35400, Q68ED2, Q14831
Entrez	108073, 81672, 2917
Expiration Date	12 months from date of receipt.