

Product Datasheet

GTF2E2 peptide (orb217380)

Description	The peptide is used to block anti-GTF2E2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GTF2E2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214024.
Purity	> 85%
Source	Synthetic
Uniprot ID	P29084, Q9D902
Entrez	68153, 2961
Expiration Date	12 months from date of receipt.