

Product Datasheet

GUCY1B3 peptide (orb217383)

Description	The peptide is used to block anti-GUCY1B3 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	GUCY1B3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214027.
Purity	> 85%
Source	Synthetic
Uniprot ID	O54865, P20595, Q02153
Entrez	25202, 54195, 2983
Expiration Date	12 months from date of receipt.