

Product Datasheet

HMOX2 peptide (orb217396)

Description	The peptide is used to block anti-Heme Oxygenase 2 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	HMOX2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214040.
Purity	> 85%
Source	Synthetic
Uniprot ID	P23711, O70252, P30519
Entrez	15369, 79239, 3163
Expiration Date	12 months from date of receipt.