

## Product Datasheet

### HNRNPA1 peptide (orb217398)

|                            |                                                                 |
|----------------------------|-----------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block anti-hnRNP A1 antibody reactivity. |
| <b>Reactivity</b>          | Human, Mouse, Porcine, Primate, Rat                             |
| <b>Tested Applications</b> | Blocking                                                        |
| <b>Target</b>              | HNRNPA1                                                         |
| <b>Form/Appearance</b>     | Lyophilized powder                                              |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                    |
| <b>Note</b>                | For research use only                                           |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb214042.        |
| <b>Purity</b>              | > 85%                                                           |
| <b>Source</b>              | Synthetic                                                       |
| <b>Uniprot ID</b>          | <b>P04256, P09651, P49312</b>                                   |
| <b>Entrez</b>              | <b>15382, 29578, 3178</b>                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                 |