

Product Datasheet

HOXA1 peptide (orb217402)

Description	The peptide is used to block anti-HOXA1 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	HOXA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214046.
Purity	> 85%
Source	Synthetic
Uniprot ID	P49639, P09022, O08656
Entrez	3198, 25607
Expiration Date	12 months from date of receipt.