

Product Datasheet

HOXB9 peptide (orb217408)

Description	The peptide is used to block anti-HOXB9 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	HOXB9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214052.
Purity	> 85%
Source	Synthetic
Uniprot ID	P17482, P20615
Entrez	3219, 15417
Expiration Date	12 months from date of receipt.