

Product Datasheet

NRAS peptide (orb217412)

Description	The peptide is used to block anti-NRAS/HRAS/KRAS antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	NRAS; HRAS; KRAS
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214056.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q61411, P20171, P01116, P01111, P01112, P08644, Q04970, P08556, P32883
Entrez	3265, 4893, 3845, 16653, 24605, 293621, 15461, 24525
Expiration Date	12 months from date of receipt.