

Product Datasheet

HRH2 peptide (orb217414)

Description	The peptide is used to block anti-Histamine H2 Receptor antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	HRH2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214058.
Purity	> 85%
Source	Synthetic
Uniprot ID	P25021, P97292, P25102
Entrez	3274, 25461, 15466
Expiration Date	12 months from date of receipt.