

Product Datasheet

HSD11B1 peptide (orb217415)

Description	The peptide is used to block anti-HSD11B1 antibody reactivity.
Reactivity	Human, Rabbit, Rat
Tested Applications	Blocking
Target	HSD11B1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214059.
Purity	> 85%
Source	Synthetic
Uniprot ID	P28845, P16232
Entrez	3290, 25116
Expiration Date	12 months from date of receipt.