

Product Datasheet

HSP90AA1 peptide (orb217421)

Description	The peptide is used to block anti-HSP90 alpha antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	HSP90AA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214065.
Purity	> 85%
Source	Synthetic
Uniprot ID	P07900, P82995, P07901
Entrez	15519, 3320, 299331
Expiration Date	12 months from date of receipt.