

Product Datasheet

HSPD1 peptide (orb217422)

Description	The peptide is used to block anti-HSP60 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	HSPD1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214066.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10809, P63038, P63039
Entrez	15510, 63868, 3329
Expiration Date	12 months from date of receipt.