

Product Datasheet

HTR1A peptide (orb217425)

Description	The peptide is used to block anti-5-HT1A antibody reactivity.
Reactivity	Human, Rabbit, Rat
Tested Applications	Blocking
Target	HTR1A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214069.
Purity	> 85%
Source	Synthetic
Uniprot ID	P19327, P08908
Entrez	3350, 24473
Expiration Date	12 months from date of receipt.