

Product Datasheet

HTR1D peptide (orb217426)

Description	The peptide is used to block anti-5-HT1D antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	HTR1D
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214070.
Purity	> 85%
Source	Synthetic
Uniprot ID	P28565, P28221, Q61224
Entrez	15552, 3352, 25323
Expiration Date	12 months from date of receipt.