

Product Datasheet

IAPP peptide (orb217431)

Description	The peptide is used to block anti-Amylin antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	IAPP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214075.
Purity	> 85%
Source	Synthetic
Uniprot ID	P12969, P10997, P12968
Entrez	24476, 3375, 15874
Expiration Date	12 months from date of receipt.