

Product Datasheet

IKBKB peptide (orb217447)

Description	The peptide is used to block anti-IKK beta
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	IKBKB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214091.
Purity	> 85%
Source	Synthetic
Uniprot ID	O88351, Q9QY78, O14920
Entrez	16150, 3551, 84351
Expiration Date	12 months from date of receipt.