

Product Datasheet

MYL9 peptide (orb217455)

Description	The peptide is used to block anti-MYL9
Reactivity	Human, Mouse, Porcine
Tested Applications	Blocking
Target	MYL9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214099.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9CQ19 , Q64122 , P24844
Entrez	296313 , 98932 , 10398
Expiration Date	12 months from date of receipt.