

Product Datasheet

IL12 p40 peptide (orb217460)

Description	The peptide is used to block anti-IL12 p40 antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	IL12B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214104.
Purity	> 85%
Source	Synthetic
Uniprot ID	P29460
Entrez	3593
Expiration Date	12 months from date of receipt.