

Product Datasheet

IRAK1 peptide (orb217468)

Description	The peptide is used to block anti-IRAK1
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	IRAK1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214112.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q62406, P51617
Entrez	16179, 3654
Expiration Date	12 months from date of receipt.