

Product Datasheet

IRS1 peptide (orb217471)

Description	The peptide is used to block anti-IRS1
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	IRS1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214115.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35569, P35570, P35568
Entrez	16367, 3667, 25467
Expiration Date	12 months from date of receipt.