

Product Datasheet

JAK2 peptide (orb217493)

Description	The peptide is used to block anti-JAK2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	JAK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214137.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q62120, Q62689, O60674
Entrez	24514, 16452, 3717
Expiration Date	12 months from date of receipt.