

Product Datasheet

JUP peptide (orb217501)

Description	The peptide is used to block anti-Gamma-catenin antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	JUP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214145.
Purity	> 85%
Source	Synthetic
Uniprot ID	P14923, Q6P0K8, Q02257
Entrez	3728, 81679, 16480
Expiration Date	12 months from date of receipt.