

Product Datasheet

KCNA5 peptide (orb217506)

Description	The peptide is used to block anti-Kv1.5 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	KCNA5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214150.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q61762, P19024, P22460
Entrez	3741, 25470, 16493
Expiration Date	12 months from date of receipt.