

Product Datasheet

KCNH1 peptide (orb217510)

Description	The peptide is used to block anti-Kv10.1 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	KCNH1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214154.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q63472, O95259, Q60603
Entrez	65198, 3756, 16510
Expiration Date	12 months from date of receipt.