

Product Datasheet

KCNJ10 peptide (orb217513)

Description	The peptide is used to block anti-Kir4.1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	KCNJ10
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214157.
Purity	> 85%
Source	Synthetic
Uniprot ID	P49655, Q9JM63, P78508
Entrez	16513, 3766, 29718
Expiration Date	12 months from date of receipt.