

Product Datasheet

c-Kit peptide (orb217518)

Description	c-Kit peptide
Reactivity	Human, Porcine, Primate
Tested Applications	Blocking
Target	KIT
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214162.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10721
Entrez	3815
Expiration Date	12 months from date of receipt.