

Product Datasheet

KRT7 peptide (orb217520)

Description	The peptide is used to block anti-Cytokeratin 7 antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	KRT7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214164.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6IG12 , Q9DCV7 , P08729
Entrez	110310 , 3855 , 300242
Expiration Date	12 months from date of receipt.