

Product Datasheet

KRT8 peptide (orb217522)

Description	The peptide is used to block anti-Cytokeratin 8 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	KRT8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214166.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q10758 , P05787 , P11679
Entrez	25626 , 3856 , 16691
Expiration Date	12 months from date of receipt.