

Product Datasheet

KRT17 peptide (orb217527)

Description	The peptide is used to block anti-Cytokeratin 17 antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	KRT17
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214171.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6IFU8, Q04695
Entrez	287702, 3872
Expiration Date	12 months from date of receipt.