

Product Datasheet

LGALS1 peptide (orb217541)

Description	The peptide is used to block anti-Galectin 1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	LGALS1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214185.
Purity	> 85%
Source	Synthetic
Uniprot ID	P11762 , P09382 , P16045
Entrez	56646 , 16852 , 3956
Expiration Date	12 months from date of receipt.