

Product Datasheet

Galectin 3 peptide (orb217542)

Description	Galectin 3 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	LGALS3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214186.
Purity	> 85%
Source	Synthetic
Uniprot ID	P16110, P08699, P17931
Entrez	83781, 3958
Expiration Date	12 months from date of receipt.