

Product Datasheet

LIMK1 peptide (orb217548)

Description	LIMK1 peptide
Reactivity	Human, Rat
Tested Applications	Blocking
Target	LIMK1; LIMK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214192.
Purity	> 85%
Source	Synthetic
Uniprot ID	P53667, P53669
Entrez	3984
Expiration Date	12 months from date of receipt.