

Product Datasheet

LMNA peptide (orb217555)

Description	The peptide is used to block anti-Lamin A/C antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	LMNA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214199.
Purity	> 85%
Source	Synthetic
Uniprot ID	P48678, P02545, P48679
Entrez	60374, 16905, 4000
Expiration Date	12 months from date of receipt.