

Product Datasheet

VWA5A peptide (orb217558)

Description	The peptide is used to block anti-VWA5A antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	VWA5A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214202.
Purity	> 85%
Source	Synthetic
Uniprot ID	O00534, Q99KC8
Entrez	4013, 67776
Expiration Date	12 months from date of receipt.