

Product Datasheet

MARCKS peptide (orb217564)

Description	The peptide is used to block anti-MARCKS
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	MARCKS
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214208.
Purity	> 85%
Source	Synthetic
Uniprot ID	P29966, P30009, P26645
Entrez	17118, 4082
Expiration Date	12 months from date of receipt.