

## Product Datasheet

### MARCKS peptide (orb217565)

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block anti-MARCKS antibody reactivity. |
| <b>Reactivity</b>          | Human, Mouse, Rat, Zebrafish                                  |
| <b>Tested Applications</b> | Blocking                                                      |
| <b>Target</b>              | MARCKS                                                        |
| <b>Form/Appearance</b>     | Lyophilized powder                                            |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                  |
| <b>Note</b>                | For research use only                                         |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb214209.      |
| <b>Purity</b>              | > 85%                                                         |
| <b>Source</b>              | Synthetic                                                     |
| <b>Uniprot ID</b>          | <b>P29966, P30009, P26645</b>                                 |
| <b>Entrez</b>              | <b>17118, 4082</b>                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                               |