

Product Datasheet

MAD2L1 peptide (orb217566)

Description	The peptide is used to block anti-MAD2 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	MAD2L1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214210.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13257, Q9Z1B5
Entrez	4085, 56150
Expiration Date	12 months from date of receipt.