

Product Datasheet

SMAD1 peptide (orb217568)

Description	The peptide is used to block anti-SMAD1 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	SMAD1; SMAD5; SMAD9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214212.
Purity	> 85%
Source	Synthetic
Uniprot ID	P70340, O54835, Q9JIW5, Q9R1V3, Q15797, P97454, O15198, Q99717, P97588
Entrez	4086, 4093, 4090, 55994, 59328, 17125, 85435, 17129
Expiration Date	12 months from date of receipt.