

Product Datasheet

SMAD2 peptide (orb217570)

Description	The peptide is used to block anti-SMAD2/3 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	SMAD2; SMAD3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214214.
Purity	> 85%
Source	Synthetic
Uniprot ID	P84025, P84022, O70436, Q15796
Entrez	4087, 25631, 29357, 4088
Expiration Date	12 months from date of receipt.