

Product Datasheet

Tau peptide (orb217581)

Description	Tau peptide
Reactivity	Human, Mouse, Rabbit, Zebrafish
Tested Applications	Blocking
Target	MAPT
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214225.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10637, P10636
Entrez	4137, 17762
Expiration Date	12 months from date of receipt.