

Product Datasheet

MIP peptide (orb217603)

Description	The peptide is used to block anti-Aquaporin 0 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	MIP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214247.
Purity	> 85%
Source	Synthetic
Uniprot ID	P30301, P51180, P09011
Entrez	17339, 4284
Expiration Date	12 months from date of receipt.