

Product Datasheet

FOXO4 peptide (orb217608)

Description	The peptide is used to block anti-FOXO4 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	FOXO4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214252.
Purity	> 85%
Source	Synthetic
Uniprot ID	P98177, Q9WVH3
Entrez	4303, 54601
Expiration Date	12 months from date of receipt.