

Product Datasheet

MNAT1 peptide (orb217624)

Description	The peptide is used to block anti-MNAT1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	MNAT1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214268.
Purity	> 85%
Source	Synthetic
Uniprot ID	P51949, P51948
Entrez	17420, 4331
Expiration Date	12 months from date of receipt.