

Product Datasheet

CD206 peptide (orb217628)

Description	CD206 antibody
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	MRC1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214272.
Purity	> 85%
Source	Synthetic
Uniprot ID	P22897, Q61830
Entrez	4360, 17533
Expiration Date	12 months from date of receipt.