

Product Datasheet

MRE11A peptide (orb217629)

Description	The peptide is used to block anti-MRE11 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	MRE11A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214273.
Purity	> 85%
Source	Synthetic
Uniprot ID	P49959 , Q61216 , Q9JIM0
Entrez	17535 , 64046 , 4361
Expiration Date	12 months from date of receipt.