

Product Datasheet

MYC peptide (orb217639)

Description	The peptide is used to block anti-c-Myc
Reactivity	Human, Mouse, Porcine, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	MYC
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214283.
Purity	> 85%
Source	Synthetic
Uniprot ID	P01106, P01108, P09416
Entrez	17869, 24577, 4609
Expiration Date	12 months from date of receipt.