

Product Datasheet

MyD88 peptide (orb217642)

Description	The peptide is used to block anti-MYD88 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat, Sheep
Tested Applications	Blocking
Target	MyD88
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214286.
Purity	> 85%
Source	Synthetic
Uniprot ID	P22366, Q6Y1S1, Q99836
Entrez	301059, 17874, 4615
Expiration Date	12 months from date of receipt.