

Product Datasheet

NFKB1 peptide (orb217661)

Description	The peptide is used to block anti-NF-kappaB p105
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	NFKB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214305.
Purity	> 85%
Source	Synthetic
Uniprot ID	P25799, P19838, Q63369
Entrez	18033, 4790
Expiration Date	12 months from date of receipt.