

Product Datasheet

NFKBIA peptide (orb217667)

Description	The peptide is used to block anti-IKB alpha antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	NFKBIA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214311.
Purity	> 85%
Source	Synthetic
Uniprot ID	P25963, Q9Z1E3
Entrez	4792, 18035
Expiration Date	12 months from date of receipt.